




Trapeze Cushion Dressing

The **LIGASANO® Trapeze Cushion Dressing** made of LIGASANO® orange is an innovative completion to the tried and tested Omega cushion. It is primarily used in intensive care, especially for ventilated patients in the prone position.

This padded dressing enables very stable positioning of the head without losing sight of the patient's face. The tube or tracheal cannula can be monitored at all times and corrected if necessary. Mouth and nose care is also possible without any problems. The trapezoidal shape and the option of height adjustment prevent the neck from being hyperextended.

Due to the increased number of prone positions for ventilated patients during the COVID pandemic, a patient-specific option for stabilisation and pressure reduction could be offered, which is successfully used in many intensive care units.






A practical feature of the LIGASANO® Trapeze Cushion Dressing is,

- that it rests firmly on the ground. This is made possible by the wide installation surface.
- that it is also ideal for patients who are overweight and have a short neck. The forehead lies as if in a bowl; nose, mouth, chin and neck are not covered.
- that it is adjustable in height, in four 1cm steps (forehead support height selectable from 9-12 cm)
- that it can also be adjusted three times in height and width in the forehead support area
- that LIGASANO® white is used as a dressing pad in the area of contact with the forehead. In this way, wounds in the forehead area can also be treated or, even better, avoided
- The use of LIGASANO® orange in combination with LIGASANO® white counteracts heat and moisture build-up in the contact area, which is far more than just an improvement in comfort.

Trapezer Cushion Dressing

Art. no.	UDI-DI	Sales unit
36060-006	111767494909	1 piece LIGASANO® orange unsterile Trapeze cushion dressing ● 5 pieces LIGASANO® white unsterile Skin contact dressing 30 x 12 x 1 cm ○



You can buy LIGASANO® products at pharmacies or medical supply stores. If an article is not on hand at the wholesaler, you can also order directly from us. We deliver of course also worldwide. We advise you the cost for shipping on request. Because this cost vary to size, weight and country, we can number the total cost only after a precise request or order.

LIGASANO®
wirkt einfach

LIGAMED® medical Produkte GmbH
Pfannenstielstr. 2 • D - 90556 Cadolzburg
Tel. 09103 / 2046 • Fax 09103 / 2796 • E-Mail: info@ligamed.de

www.ligasano.com